

ARMYGIL

emergency and field army hospitals
X-ray mobile unit



GILARDONI

Mobile radiological unit specifically designed for use in Field Army Hospitals and in any emergency situation.

Armygil is dedicated for the execution of a complete radiographic and fluoroscopic examination.

The system provides wide vertical adjustment (focus at 2 meters from the ground), for radiographic imaging needing an high focus-film distance, also on stretcher-case patients.

Large diameter wheels allow use of the equipment also on uneven grounds.

Armygil is completely disassemblable without the need of tools and can be packed in cases, easily transportable on pallets.

A foldable patient stretcher, provided with horizontal and vertical cassette holder, allows projections both with patient in orthostatic and troco position.

Technical Features:

X-ray generator:	high frequency, 40 kHz
X-ray tube:	rotating anode
Focus sizes:	0,3 - 0,6 mm
kV range:	40 / 120 kV
Max mA (radiography mode):	65 mA
Max mA (fluoroscopy mode):	4 mA
Image Intensifier:	7/5" (18/12 cm)
Video system:	CCD - 625 lines
Memory (basic configuration):	1+48 images
Monitors:	2, size 17"
Cassette holder:	24x30 cm



Reserving the right to improve and modify

Acknowledged as "HIGHLY QUALIFIED LABORATORY" with decree D.M. 9-10-1985 - L.46/82 art.4

Head office and facilities: Via **Arturo Gilardoni**, 1 - 23826 Mandello del Lario (LC) - Italy
tel. (+39) 0341.705.111 - fax (+39) 0341.735.046
e-mail: gx@gilardoni.it - www.gilardoni.it

Export department: tel. (+39) 0341.705.225 - 0341.705.241

